SINGLE-COLOR ANTIBODIES

Anti-ZAP-70-PE

PN B10248

AntibodyReagentSpecificity:HumanLine:Clone:SBZAPFormat:Isotype:IgG1 mouseSize:

Line: IOTest Format: Liquid Size: 1 mL Status: ASR

Antibody characteristics

The monoclonal antibody SBZAP was generated against a KLH-peptide sequence corresponding to the human ZAP-70 amino acid residues 280-309.


Antigen characteristics

ZAP-70, a Syk family protein tyrosine kinase expressed in T and NK cells, plays a critical role in mediating T cell activation in response to T cell receptor (TCR) engagement. Following TCR engagement, ZAP-70 is rapidly phosphorylated on several tyrosine residues, presumably by two mechanisms: autophosphorylation and transposphorylation by the Src family tyrosine kinase, Lck. Tyrosine phosphorylation of ZAP-70 correlates with its increased kinase activity and downstream signalling events.

Formats available Size Line Status Part # PE 1 mL IOTest ASR B10248

ASR: Analyte Specific Reagent. Analytical and performance characteristics are not established.

Example of results



NEW PRODUCT

Permeabilized whole blood sample

References

- Iwashima M., Irving B.A., van Oers N.S., Chan A.C. and Weiss A. 1994. "Sequential interactions of the TCR with two distinct cytoplasmic tyrosine kinases." Science, 263, 1136-1139.
- Neumeister E.N., Zhu Y., Richard S., Tershorst C. Chan A.C. and Shaw A.S. 1995. "Binding of ZAP-70 to phosphorylated T-cell receptor zeta and eta enhances its autophosphorylation and generates specific binding sites for SH2 domaincontaining proteins." Mol. Cell. Biol., 15, 3171-3178.
- Chan A.C., Dalton M., Johnson R., Kong G.H., Wang T., Thoma R. and Kurosaki T., 1995. "Activation of ZAP-70 kinase activity by phosphorylation of tyrosine 493 is required for lymphocyte antigen receptor function." EMBO J., 14, 2499-2508.
- Wange R.L., Isakov N., Burke T.R., Otaka A., Roller P.P., Watts J.D., Aebersold R. and Samelson L.E. 1995. "F2(Pmp)2-TAM zeta 3, a novel competitive inhibitor of the binding of ZAP-70 to the T cell antigen receptor, blocks early T cell signaling." J. Biol. Chem. 270, 944.

- Van Oers, N.S., and A. Weiss. 1995. "The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes." Semin. Immunol. 7 (4), 227-236.
- CHU D.H., Morita C.T. and Weiss A. 1998. "The Syk family of protein tyrosine kinases in T-cell activation and development." Immunol. Rev.165, 167-180.
- William B.L., Irvin B.J., Sutor S.L., Chini C.C., Yacyshyn E., Bubeck W.J., Dalton M., Chan A.C. and Abraham R.T. 1999.
 "Phosphorylation of Tyr319 in ZAP-70 is required for T-cell antigen receptordependant phospholipase C-gamma1 and Ras activation." EMBO J. 18, 1832-1844.
- Di Bartolo V., Mége D. Germain V., Pelosi M., Dufour E., Michel F., Magistrelli G., Isacchi A. and Acuto O. 1999. "Tyrosine 319, a newly indentified phosphorylation site of ZAP-70, plays a critical role in T cell antigen receptor signaling." J. Biol. Chem., 274, 6285-6294.

BECKMAN COULTER Life Sciences